Surgical Gown CE/FDA/TEST REPORT Approved









Certificate

No. Q6 083528 0010 Rev. 00

Holder of Certificate: Henan JoinKona Medical

Products Stock Co., Ltd.

Xinxing Road

The South of Industry District

LuShan County

467300 PingDingShan, Henan Province PEOPLE'S REPUBLIC OF CHINA

Facility(ies): Henan JoinKona Medical Products Stock Co., Ltd.

Xinxing Road, The South of Industry District, LuShan County, 467300 PingDingShan, Henan Province, PEOPLE'S REPUBLIC

OF CHINA

Certification Mark:



Production and Distribution of Surgical Set, Scope of Certificate:

Surgical Gown, Surgical Drape, Tube Cover, Liquid Collection Pouch, Protective Gown,

Warm Blanket and Surgical Hood

Applied Standard(s): EN ISO 13485:2016

Medical devices - Quality management systems -

Requirements for regulatory purposes

(ISO 13485:2016) DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system (excluding subclause 7.3), which meets the requirements of the listed standard(s). See also notes overleaf.

SH1978307 Report No .:

2019-04-01 Valid from: Valid until: 2022-03-31

Date, 2019-04-01

1. Pumil



EC Certificate

Production Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex V (Devices in class I in sterile conditions, sterilised systems or procedure packs)

No. G2S 18 01 83528 008

Manufacturer: Henan JoinKona Medical

Products Stock Co., Ltd.

Xinxing Road

The South of Industry District

LuShan County

467300 PingDingShan, Henan Province PEOPLE'S REPUBLIC OF CHINA

EC-Representative: Shanghai International Holding

Corp. GmbH (Europe)

Eiffestraße 80 20537 Hamburg GERMANY

Product

Category(ies):

Surgical Set, Surgical Gown, Surgical Drape, Tube Cover, Liquid Collection Pouch,

Warm Blanket and Surgical Hood

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for manufacture in accordance with MDD Annex V. This quality assurance system covers those aspects of manufacture concerned with securing and maintaining sterile conditions of the respective devices / device categories and conforms to the requirements of this Directive. It is subject to periodical surveillance. See also notes overleaf.

Report No.: SH18783EXT01

Valid from: 2018-04-23 Valid until: 2023-04-22

Date, 2018-02-15

Stefan Preiß







Fiscal Year 2020 CERTIFICATION OF REGISTRATION

This certifies that:

HENAN JOINKONA MEDICAL PRODUCTS STOCK CO.,LTD

Xinxing Road, The South of Industry District, Lushan County, Pingdingshan, HENAN, 467300, CHINA

has completed the FDA Establishment Registration and Device Listing with the US Food & Drug Administration, through

Shenzhen CTB Testing Technology Co., Ltd.

Owner/Operator Number: 10062713

Device Listing#: See annex

CTB will confirm that such registration remains effective upon request and presentation of this certificate until the end of the calendar year stated above, unless said registration is terminated after issuance of this certificate. CTB makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. CTB assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding." The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration, CTB is not affiliated with the U.S. Food and Drug Administration.





Issued: March 4, 2020

Expiration Date: December 31, 2020



Fiscal Year 2020 CERTIFICATION OF REGISTRATION

Annex to Device Listing# for Owner/Operator Number: 10062713

Listing No.	Code	Device Name	Proprietary Names	Activities
47		INSTRUMENT,	FEMALE SS	
D373299	MDM	MANUAL,	3.5MM PELLET	
D373299	MIDIVI	SURGICAL,	INSERTION	
		GENERAL USE	TRAY	
D373300	KDD	Kit, surgical instrument.	Surgical Set	
D373300	KDD	disposable	Surgical Set	
D373303	FYE	DRESS, SURGICAL	Surgical Gown	Manufacturer Contract
D373304	KME	BEDDING, DISPOSABLE, MEDICAL	Surgical Drape	Manufacturer Contract Sterilizer Repackager/Relabeler
D373305	FMW	COVER, MATTRESS (MEDICAL PURPOSES)	Tube Cover	Remanufacturer Foreign Exporter
D373306	KET	FILTERS, CELL COLLECTION, TISSUE PROCESSING	Liquid Collection Pouch	
D373307	FXZ	HELMET, SURGICAL	Warm Blanket and Surgical Hood	

NOT END OF THE ANNEX



Fiscal Year 2020 CERTIFICATION OF REGISTRATION

Annex to Device Listing# for Owner/Operator Number: 10062713

Listing No.	Code	Device Name	Proprietary Names	Activities
D373308	OEA	Non-surgical isolation gown	Isolation Gown	Manufacturer Contract
D373309	KPY	Shield, protective, personnel	Protective Coverall	Manufacturer Contract Sterilizer Repackager/Relabeler Remanufacturer Foreign Exporter
D373310	IMD	PACK, HOT OR COLD, DISPOSABLE	Dressing Pack	

END OF THE ANNEX

Issued: March 4, 2020 Expiration Date: December 31, 2020





TEST REPORT

Report No.:

SHT16110079-01E

Date:

2016-12-05

Page 1of 2

The following information of testing sample

Sumple Name	Disposable Surgical Gown	Sample No.	116110079-01
Sample Type	1	Sample Qty.	4pcs
Applicant	Henan Joi	inkona Medical Produ	cts Stock Co.,Ltd
Address	Xinxing Road, The South of Ir	ndustry District, LuSh Province, Chin	an County, PingDingShan City, a
Sample Received Date	2016-10-28 Testing period 2016-		2016-10-28 to 2016-12-01
Test Items	Water	resistance, Hydrostati	ic pressure test
Test Method	AATCC 42-2013 AATCC 127-2013		

*****For more detailed information, please refer to next page(s)*****

Redacted by:

(Yufang Wang)

(Tao Wu)

Approved by:

(Authorized signatory: Fei Xia)





TEST REPORT

Report No.:

SH1116 (10079-011)

Date:

2016-12-05

Page 2of 2

Test Result(s):

I.Water resistance

Test Method: AATCC 42-2013

Determinate Standard: According to applicant's requirement

Test Item	Unit	Test Result	Limited Value	Determination	
Water resistance	g	0.4	≤1.0	Pass	

2.Hydrostatic pressure test

Test Method: AATCC 127-2013

Determinate Standard: According to applicant's requirement

Test Item	Unit	Test Result	Limited Value	Determination
Hydrostatic pressure test	cmH ₂ O	65.2	≥50	Pass

Sample Photo



The End

Remark: The samples provided by the customer, the report test results are only responsible for samples.

The test data is only used for reference, not as a social justice data.

Disposable Surgical Gown

1. Description

Disposable Surgical Gowns are critical items of protective clothing suitable for Operation Room, medical clinics, hospital ward, inspection rooms, laboratories, ICU and CDC sites for important isolation of virus damage.

At Joinkona Medical Company, there is a wide selection of disposable surgical gowns made of SMMS that will certainly serve to protect medical staffs in situations where there are exposure concerns.

2. Advantages to wearing disposable Surgical Gown

Raw Material:

Weight: 35g/M² or 45g/M² SMMS Non-woven Fabrics.

SMMS Performance: Fluid Proof and Anti-Static or Fluid Proof and Anti-Static and Anti-Alcohol and Anti-Blood(AS&AR)

In the healthcare profession, disposable surgical gown play a crucial role in asepsis by reducing the transfer of bacteria from the skin of the medical staff to the air and protect medical staff in operation room. Besides, It will also protect staff members from blood, urine, saline, or other chemicals and bodily fluids during surgical procedures.

Comfortable, Lightweight and Breathable

They are breathable allowing the wearer to be insulated and not overheat

3. Specifications and Features

Style

One piece of Disposable Surgical Gown

Features

- Single use
- Secure protection (ultrasonic technology)
- Anti-fluid, Anti-alcohol, Anti-blood, Anti-static
- Durable
- Comfortable, Lightweight and Breathable
- > Tear-resistant
- Flame retardant






- > CE
- > FDA
- > ISO 13485
- > EN13795

Standard Surgical Gown

Reinforced Surgical Gown



Surgical Gown Size					Unit:cm		
Size	s	M	L	XL	XXL	XXXL	Tolerance
A	110	115	127	135	140	145	±2
B	140	145	150	150	155	160	±2

4. Workmanship Details

- 1, Ultrasonic technology in arm and sleeve.
- 2, Sleeve with cuff.
- 3, Four belts with belt card .
- 4, comfortable for collar with soft white spunlace.

5. Products Package:

1PC/Sterile Pouch

35gsm SMMS Standard Surgical Gown and Reinforced Surgical Gown

Products Name	Size	Carton Size	Qty/Carton (pcs)
Standard Surgical Gown	S	50x40x45cm	50
Standard Surgical Gown	M	50x40x45cm	50
Standard Surgical Gown	L	50x40x45cm	50
Standard Surgical Gown	XL	50x40x45cm	50
Standard Surgical Gown	2XL	50x40x45cm	50
Standard Surgical Gown	3XL	50x40x45cm	50
Reinforced Surgical Gown	S	50x40x45cm	50
Reinforced Surgical Gown	M	50x40x45cm	50
Reinforced Surgical Gown	L	50x40x45cm	50
Reinforced Surgical Gown	XL	50x40x45cm	50
Reinforced Surgical Gown	2,/L	50x40x45cm	50
Reinforced Surgical Gown	3XL	50x40x45cm	50

45gsm SMMS Standard Surgical Cown and Reinforced Surgical Gown

Products Name	Size	Carton Size	Qty/Carton (pcs)
Standard Surgical Gown	S	50x40x45cm	40
Standard Surgical Gown	M	50x40x45cm	40
Standard Surg cal Gown	L	50x40x45cm	40
Standard Surgical Gown	XL	50x40x45cm	40
Standard Surgical Gown	2XL	50x40x45cm	40
Standard Surgical Gown	3XL	50x40x45cm	40
Reinforced Surgical Gown	S	50x40x45cm	40
Reinforced Surgical Gown	М	50x40x45cm	40
Reinforced Surgical Gown	L	50x40x45cm	40
Reinforced Surgical Gown	XL	50x40x45cm	40
Reinforced Surgical Gown	2XL	50x40x45cm	40
Reinforced Surgical Gown	3XL	50x40x45cm	40











